



Technical Specification

ISO/TS 20224-10

Molecular biomarker analysis — Detection of animal-derived materials in foodstuffs and feedstuffs by real-time PCR —

Part 10: Duck DNA detection method

*Analyse de biomarqueurs moléculaires — Détection de matériaux
d'origine animale dans les denrées alimentaires et les aliments
pour animaux par PCR en temps réel —*

Partie 10: Méthode de détection de l'ADN de canard

**First edition
2024-02**

STANDARDSISO.COM : Click to view the full PDF of ISO/TS 20224-10:2024



COPYRIGHT PROTECTED DOCUMENT

© ISO 2024

All rights reserved. Unless otherwise specified, or required in the context of its implementation, no part of this publication may be reproduced or utilized otherwise in any form or by any means, electronic or mechanical, including photocopying, or posting on the internet or an intranet, without prior written permission. Permission can be requested from either ISO at the address below or ISO's member body in the country of the requester.

ISO copyright office
CP 401 • Ch. de Blandonnet 8
CH-1214 Vernier, Geneva
Phone: +41 22 749 01 11
Email: copyright@iso.org
Website: www.iso.org

Published in Switzerland

Contents

Page

Foreword	iv
Introduction	v
1 Scope	1
2 Normative references	1
3 Terms and definitions	1
4 Scientific basis	2
5 Reagents and materials	2
5.1 General	2
5.2 PCR reagents	2
6 Apparatus	3
7 Procedure	3
7.1 Preparation of the test portion/sample	3
7.2 Preparation of DNA extracts	3
7.3 PCR setup	3
7.3.1 Reaction mixes	3
7.3.2 PCR controls	4
7.3.3 Real-time PCR thermocycler plate set-up	4
7.4 Temperature-time programme	4
8 Accept/reject criteria	5
8.1 General	5
8.2 Identification	5
9 Validation status and performance criteria	5
9.1 General	5
9.2 Robustness	5
9.3 Reproducibility	6
9.4 Sensitivity	6
9.5 Specificity	9
10 Test report	11
Annex A (informative) BlastN +2.12.0 results for query of GenBank RefSeq genome (refseq_genomes) and whole-genome shotgun contigs (wgs)	12
Annex B (informative) Members of the <i>Anatidae</i> family and its family tree established with available public genomic sequences	16
Bibliography	23

Foreword

ISO (the International Organization for Standardization) is a worldwide federation of national standards bodies (ISO member bodies). The work of preparing International Standards is normally carried out through ISO technical committees. Each member body interested in a subject for which a technical committee has been established has the right to be represented on that committee. International organizations, governmental and non-governmental, in liaison with ISO, also take part in the work. ISO collaborates closely with the International Electrotechnical Commission (IEC) on all matters of electrotechnical standardization.

The procedures used to develop this document and those intended for its further maintenance are described in the ISO/IEC Directives, Part 1. In particular, the different approval criteria needed for the different types of ISO document should be noted. This document was drafted in accordance with the editorial rules of the ISO/IEC Directives, Part 2 (see www.iso.org/directives).

ISO draws attention to the possibility that the implementation of this document may involve the use of (a) patent(s). ISO takes no position concerning the evidence, validity or applicability of any claimed patent rights in respect thereof. As of the date of publication of this document, ISO had not received notice of (a) patent(s) which may be required to implement this document. However, implementers are cautioned that this may not represent the latest information, which may be obtained from the patent database available at www.iso.org/patents. ISO shall not be held responsible for identifying any or all such patent rights.

Any trade name used in this document is information given for the convenience of users and does not constitute an endorsement.

For an explanation of the voluntary nature of standards, the meaning of ISO specific terms and expressions related to conformity assessment, as well as information about ISO's adherence to the World Trade Organization (WTO) principles in the Technical Barriers to Trade (TBT), see www.iso.org/iso/foreword.html.

This document was prepared by Technical Committee ISO/TC 34, *Food products*, Subcommittee SC 16, *Horizontal methods for molecular biomarker analysis*.

A list of all parts in the ISO 20224 series can be found on the ISO website.

Any feedback or questions on this document should be directed to the user's national standards body. A complete listing of these bodies can be found at www.iso.org/members.html.

Introduction

Fraudulent adulteration of meat in food and feed threatens both public safety and commerce. Adulteration can affect those adhering to ethnological dietary rules, economic development and social stability. This document provides a real-time polymerase chain reaction (real-time PCR) analytical method for the identification of meat animal species from nucleic acid present in the ingredients of food and feed.

Animal-derived biological materials in food and feed are detected and identified in the laboratory with the following successive (or simultaneous) steps: preparation of the test portion/sample, nucleic acid extraction and purification, PCR amplification and interpretation of results. This document provides guidance for PCR amplification and interpretation of results, specific to mallard duck (*Anas platyrhynchos*) and spot-billed duck (*Anas zonorhyncha*) DNA detection. Cross detection of white-winged duck (*Asarcornis scutulata*), tufted duck (*Aythya fuligula*), muscovy duck (*Cairina moschata*) and Mandarin duck (*Aix galericulata*) is observed.

The ISO 20224 series consists of technical specifications that describe specific applications. New species DNA detection methods established in the future will be added as independent parts.

STANDARDSISO.COM : Click to view the full PDF of ISO/TS 20224-10:2024

STANDARDSISO.COM : Click to view the full PDF of ISO/TS 20224-10:2024

Molecular biomarker analysis — Detection of animal-derived materials in foodstuffs and feedstuffs by real-time PCR —

Part 10: Duck DNA detection method

1 Scope

This document specifies a real-time polymerase chain reaction (real-time PCR) method for the qualitative detection of duck-specific DNA derived from food and feed. It requires the extraction of an adequate amount of PCR amplifiable DNA from the relevant matrix and can be applied to the detection of duck material derived from mallard duck (*Anas platyrhynchos*) and spot-billed duck (*Anas zonorhyncha*). Cross detection of white-winged duck (*Asarcornis scutulata*), tufted duck (*Aythya fuligula*), muscovy duck (*Cairina moschata*) and Mandarin duck (*Aix galericulata*) is observed. Mallard duck and muscovy duck are domesticated poultry for food, while spot-billed duck, white-winged duck, tufted duck and Mandarin duck are wild avian.

The target sequence is a partial fragment of the *Anas platyrhynchos* breed pekin duck isolate CAU_Pekin_2.0 Chr1, whole genome shotgun sequence (i.e. GenBank accession number JACEUM010000001.1),^[1] which is present as a single copy per haploid genome. The provided PCR assay for this target has an absolute limit of detection of five copies per reaction, with ≥ 95 % confidence at this concentration ($\text{LOD}_{95\%}$).

2 Normative references

The following documents are referred to in the text in such a way that some or all of their content constitutes requirements of this document. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

ISO 16577, *Molecular biomarker analysis — Vocabulary for molecular biomarker analytical methods in agriculture and food production*

ISO 20813, *Molecular biomarker analysis — Methods of analysis for the detection and identification of animal species in foods and food products (nucleic acid-based methods) — General requirements and definitions*

ISO 21571, *Foodstuffs — Methods of analysis for the detection of genetically modified organisms and derived products — Nucleic acid extraction*

ISO 24276, *Foodstuffs — Methods of analysis for the detection of genetically modified organisms and derived products — General requirements and definitions*

3 Terms and definitions

For the purposes of this document, the terms and definitions given in ISO 16577 apply.

ISO and IEC maintain terminology databases for use in standardization at the following addresses:

- ISO Online browsing platform: available at <https://www.iso.org/obp>
- IEC Electropedia: available at <https://www.electropedia.org/>

4 Scientific basis

DNA is extracted from the test portion by applying a suitable method (see ISO 21571:2005, A.1). The DNA analysis consists of two parts:

- verification of the quality and amplifiability of the extracted DNA using a PCR assay specific for eukaryotes (i.e. 18S rRNA gene) or mammals and poultry (i.e. myostatin gene);
- detection of the duck species-specific DNA sequence of the single-copy *Anas platyrhynchos* breed pekin duck isolate CAU_Pekin_2.0 Chr1, whole genome shotgun sequence (i.e. GenBank accession number JACEUM010000001.1) in a real-time PCR.

NOTE The copy number of the eukaryotic ribosomal 18S RNA (18S rRNA) gene in a cell varies from several hundred to several thousand, while the specific target sequence in the duck genome and myostatin gene in mammals and poultry genome are single copy. The copy number of the specific target sequence in *Anas platyrhynchos* genome was confirmed by bioinformatics analysis at the whole genome scale (see [Annex A](#)) and digital PCR for absolute quantification.

5 Reagents and materials

5.1 General

For this document, only reagents and water of recognized analytical grade, appropriate for molecular biology, shall be used. Unless stated otherwise, solutions should be prepared by dissolving the corresponding reagents in water followed by autoclave sterilization. For all operations in which gloves are used, gloves shall be powder free. The use of aerosol protected pipette tips (protection against cross-contamination) is recommended.

5.2 PCR reagents

5.2.1 PCR master mix.

In general, real-time PCR master mix contains thermostable DNA polymerase, dNTPs, MgCl₂, KCl, and buffer as a dilutable concentrated mixture, that is ready to use.

NOTE The commercial real-time PCR master mix can be used.

5.2.2 Oligonucleotides.

The quality of the oligonucleotides shall be sufficient for their use as primers and probes. See [Table 1](#).

Table 1 — Oligonucleotides

Name	DNA sequence of the oligonucleotide	Final concentration in PCR
Specific sequence in <i>Anas platyrhynchos</i> breed pekin duck isolate CAU_Pekin_2.0 Chr1, whole genome shotgun sequence (i.e. GenBank accession number JACEUM010000001.1) ^a		
Duck-105bp-F	5'-TCTTCACAAGCAGGGTCTAATGG-3'	400 nmol/l
Duck-105bp-R	5'-CTTGGCAGAAGGTCCAGAGG-3'	400 nmol/l
Duck-105bp-P	5'- [FAM]-AGGCACAGCACGCATCTCACCACA-[TAMRA] ^b -3'	200 nmol/l
^a PCR product = 200 064 152 - TCTTCACAAG CAGGGTCTAA TGGAAGACTT GCTGGCCTGC TCTCTACTGG TGATGATGTG GTGAGATGCG TGCTGTGCCT CTCTTCCTCT GGACCTTCTG CCAAG - 200 064 048 - JACEUM010000001.1		
^b FAM: 6-carboxyfluorescein, TAMRA: 6-carboxytetramethylrhodamine.		

Duck-105bp-F is base pairs 200 064 152 – 200 064 130, Duck-105bp-R is base pairs 200 064 048 – 200 064 067 and Duck-105bp-P is 200 064 073 – 200 064 096 of JACEUM010000001.1, *Anas platyrhynchos* breed pekin duck isolate CAU_Pekin_2.0 Chr1, whole genome shotgun sequence. Equivalent reporter dyes and/or quencher dyes can be used if they yield the same or better results.

6 Apparatus

Requirements concerning apparatus and materials shall follow ISO 20813. In addition to the usual laboratory equipment, the following equipment is required.

6.1 Real-time thermocycler instrument.

A device that amplifies DNA *in vitro* and performs the temperature-time cycles is needed for PCR. Additionally, the device shall be capable of exciting fluorescence molecules at specific wavelengths and detecting sufficient emitted fluorescent light of the fluorophore used to perform TaqMan format assays.

7 Procedure

7.1 Preparation of the test portion/sample

The test sample used for DNA extraction shall be representative of the laboratory sample and homogeneous, e.g. by grinding or homogenizing the laboratory sample to a fine mixture. Test portion/sample preparation shall follow the general requirements and specific methods described in ISO 21571 and ISO 20813.

7.2 Preparation of DNA extracts

The extraction/purification and quantification of DNA from the test portion shall follow the general requirements and methods provided in ISO 21571. DNA extraction methods described in ISO 21571:2005, Annex A, are recommended.

7.3 PCR setup

7.3.1 Reaction mixes

The method is for a total volume of 25 µl per PCR. The reaction setup is given in [Table 2](#). Reagents shall be completely thawed at room temperature. Each reagent shall be carefully mixed and briefly centrifuged immediately before pipetting. A PCR reagent mixture is prepared to contain all components except for the sample DNA. The required total amount of the PCR reagent mixture prepared depends on the number of reactions to be performed, including at least one additional reaction as a pipetting reserve. The number of sample and control replicates shall follow ISO 20813. Set up the PCR tests as follows:

- mix the PCR reagent mixture, centrifuge briefly and pipette 20 µl into each reaction vial;
- add 5 µl of each sample DNA or positive DNA target control or extraction blank control or water to the respective reaction vials;
- mix and centrifuge briefly.

Table 2 — Reaction setup for the amplification

Total reaction volume	25 µl
Sample DNA ^a or controls	5 µl
2 × PCR master mix ^b	12,5 µl
Primer Duck-105bp-F, c = 10 µmol/l and Duck-105bp-R, c = 10 µmol/l	1,0 µl for each
Probe Duck-105bp-P, c = 10 µmol/l	0,5 µl
Water	to 25 µl
^a The amount of DNA in one reaction can be up to 200 ng, but the recommendation amount is less than 200 ng per reaction.	
^b In the collaborative trial, a ready-to-use optimized 2 × PCR master mix containing all of the components, excluding the template and primers, was used. The 2 × PCR master mix contains thermostable DNA polymerase, a blend of dNTPs with dUTP and uracil-UDG to minimize carry-over PCR contamination, and a passive internal reference based on ROX dye. Equivalent products can be used if they yield the same or better results. If necessary, the amounts of the reagents and the temperature-time programme can be adapted.	

7.3.2 PCR controls

7.3.2.1 General

PCR controls shall be as described in ISO 24276 and ISO 20813.

7.3.2.2 Inhibition control (reference gene assay)

A reference control gene (i.e. 18S rRNA gene for eukaryotes, myostatin gene for mammals and poultry) PCR assay using sample DNAs shall be performed to test nucleic acid amplifiability and provide control to exclude false-negative results.

7.3.3 Real-time PCR thermocycler plate set-up

Transfer the setup reaction vials to the thermocycler. The vials should be arranged to avoid any possible edge temperature variations associated with a particular real-time thermocycler instrument. Start the temperature-time programme.

7.4 Temperature-time programme

The temperature-time programme as outlined in [Table 3](#) was used in the validation study. The use of different reaction conditions and real-time PCR cycles shall be verified. The time for initial denaturation depends on the master mix used.

Table 3 — Temperature-time programme

Step	Parameter		Temperature	Time	Fluorescence measurement	Cycles
1	UNG activation ^a		50 °C	2 min	no	1
2	Initial denaturation		95 °C	10 min	no	1
3	Amplification	Denaturation	95 °C	15 s	no	45
		Annealing and elongation	60 °C	60 s	yes	
^a UNG(Uracil-N-Glycosylase) activation is mandatory if UDG-glycosylase is included in mastermix and optional if UDG-glycosylase is not included in mastermix.						

8 Accept/reject criteria

8.1 General

A corresponding real-time PCR-instrument-specific data analysis programme shall be used for the identification of PCR products. The amplification results can be expressed differently, depending on the instrument used. In the absence of detectable PCR products (e.g. negative controls), the result shall be expressed as “undetermined”, “no amplification” or the maximum number of reaction cycles performed. If amplification of the DNA target sequence in a sample (e.g. positive controls) occurred, a sigmoid-shaped amplification curve shall be observed. The cycle number at the crossing point of the amplification curve and the fluorescence threshold shall be calculated [cycle threshold (C_t) or cycle quantification (C_q)].

If, due to atypical fluorescence measurement data, the automatic interpretation does not provide a meaningful result, it can be necessary to set the baseline and the threshold manually prior to interpreting the data. In such a case, the device-specific instructions provided with the interpretation software shall be followed.

8.2 Identification

The target sequence is considered as detected if:

- duck-specific primers Duck-105bp-F and Duck-105bp-R and the probe Duck-105bp-P, produce a sigmoid-shaped amplification curve and a C_t value or C_q value $\leq \text{LOD}_{95\%}$ can be calculated;
- PCR control reactions with no added DNA (PCR reagent control, extraction blank control) produce no amplification;
- the amplification controls (positive DNA target control, PCR inhibition control) produce the expected amplification and C_t values (or C_q values).

Trace detections are defined as PCRs with C_t values later than that defined at the target $\text{LOD}_{95\%}$. In the event of a trace detection or contradictory positive/negative results from different extracts of the same sample, then the sample shall be retested. At least two new extracts shall be prepared from the homogenized laboratory sample. A minimum of 20 PCR replicates shall be conducted across the new extracts (e.g. ten PCR repeats for two extracted DNA, seven PCR repeats for three extracted DNA). The target sequence shall be considered as “detected” if $\geq 95\%$ of the new extract PCR results show a positive detection. The target sequence shall be considered as “not detected” if $< 95\%$ of the new extract PCR results show a positive detection.

9 Validation status and performance criteria

9.1 General

Validation followed a two-part process:

- a) in-house validation;
- b) collaborative trial validation.

9.2 Robustness

The robustness of the method was confirmed for the collaborative trial by changing the reaction conditions for the following factors:

- a) real-time PCR instruments (e.g. ABI 7500, BioRad CFX96, ABI 7900 HT Fast, Eppendorf Realplex 4¹⁾);

1) These are examples of a suitable products available commercially. This information is given for the convenience of users of this document and does not constitute an endorsement by ISO of these products.

- b) reaction volume: 19 µl or 21 µl PCR reagent mixture plus 5 µl sample DNA;
- c) annealing temperature: 59 °C and 61 °C;
- d) primer or probe concentration: both reduced by 30 %.

For each factor tested, the PCRs were analysed in triplicate, each with 20 copies of the target sequence and with 100 copies of the non-target sequence as negative controls. Method performance remained satisfactory for both samples and negative controls under the changed conditions for each changed factor.

9.3 Reproducibility

The reproducibility of the method was verified in a collaborative trial with 12 participants, organized by the Technical Center for Animal, Plant and Food Inspection and Quarantine, Shanghai Customs in accordance with the IUPAC protocol^[2] and the BVL guidelines.^[3] Participants received 12 DNA samples for the evaluation of false-positive and false-negative rates. All samples were labelled with randomized coding numbers and consisted of six replicate samples. The 12 DNA samples were:

- six vials of duck DNA solution, 10 copies/µl;
- six vials of bovine DNA solution, 20 copies/µl.

The copy numbers were determined using the real-time PCR of this method and serial dilutions of plasmid DNA containing the target sequence. The concentration of the plasmid DNA (copies/µl) was measured by digital PCR.

Participants received a PCR master mix and the oligonucleotides (primers and probes) from the collaborative test organizer to conduct the PCR experiments.

Duck and bovine genomic DNA were extracted from duck meat and bovine meat, respectively, and then adjusted with 0,2 × TE buffer to a nominal concentration of 10 copies/µl for duck DNA and 20 copies/µl for bovine DNA, respectively.

The collaborative trial was designed to determine false-positive and false-negative rates. Each DNA sample was tested by the participants in a single PCR test with 5 µl of the respective DNA solution, using the procedure and the conditions given in [Tables 2](#) and [3](#). The results of the collaborative trial are listed in [Table 4](#).


Table 4 — Results of the collaborative trial

Number of laboratories	12
Number of laboratories submitting results	12
Number of samples per laboratory	12
Number of accepted results	144
Number of accepted samples containing duck material	72
Number of accepted samples not containing duck material	72
False-positive results	0
False-positive results (in %)	0
False-negative results	0
False-negative results (in %)	0

9.4 Sensitivity

The absolute limit of detection (LOD_{95 %}) for the method is five DNA copies. The collaborative trial of the duck detection method was carried out at the same time as collaborative trials for the goose, turkey, and pigeon detection methods. Goose (*Anser anser domesticus*), duck (*Anas platyrhynchos*), turkey (*Meleagris gallopavo*) and rock pigeon (*Columba livia*) target DNA sequences were synthesized and cloned into the pUC57 vector (2 710 bp in length, GenBank/EMBL accession number Y14837). This constructed plasmid

pUC57-gdtp (3 167 bp in length) was sequenced to ensure that only one copy of the goose, duck, turkey and rock pigeon target DNA sequence was inserted (see [Figure 1](#)). No deletion or insertion mutations were found in the inserted sequences (see [Figure 2](#)). The target sequences of corresponding PCR methods are indicated.



Key

- 1 nt 1-121 = goose amplicon (121 bp)
- 2 nt 122~226 = duck amplicon (105 bp)
- 3 nt 227~344 = turkey amplicon (118 bp)
- 4 nt 345~457 = pigeon amplicon (113 bp)
- 5 M13 reverse promoter
- 6 ColE1 origin of replication
- 7 β -lactamase gene (ampicillin resistance gene)
- 8 ampicillin resistance gene promoter
- 9 M13 forward promoter

Figure 1 — Map of the multi-target DNA plasmid

1	<u>ACGAGGATAG</u>	<u>GTTGTGACAG</u>	<u>CTGACTCTGT</u>	<u>TCAGCCTTGC</u>	<u>GAAGACCTTA</u>
51	<u>TGCTGTCTAC</u>	<u>AATTACCTAA</u>	<u>TTGGAGGATA</u>	<u>TAGAATTATA</u>	<u>GAATCATATA</u>
101	<u>GAGAAGACGA</u>	<u>CACAGAGATT</u>	<u>CTCTTCACAA</u>	<u>GCAGGGTCTA</u>	<u>ATGGAAGACT</u>
151	<u>TGCTGGCCTG</u>	<u>CTCTCTACTG</u>	<u>GTGATGATGT</u>	<u>GGTGAGATGC</u>	<u>GTGCTGTGCC</u>
201	<u>TCTCTTCCTC</u>	<u>TGGACCTTCT</u>	<u>GCCAAGTGAA</u>	<u>CAAAATCCACT</u>	<u>TCCCTTTAAC</u>
251	<u>CTCAGGAACA</u>	<u>TCCAGCATAT</u>	<u>TGGTAAACAG</u>	<u>CGGGATGTGG</u>	<u>GGGTGTGGCT</u>
301	<u>GCGGCTCGTC</u>	<u>ATCACCTGCA</u>	<u>GCTCACTTGT</u>	<u>GCAGCAGCTA</u>	<u>ATGAGCAGTT</u>
351	<u>GTTTAGTCCT</u>	<u>CCTGTAACAC</u>	<u>GGACTCCTAA</u>	<u>GAGCACTTCT</u>	<u>CAGCCTGGCT</u>
401	<u>TTGTTTTCGT</u>	<u>CACACTGTGT</u>	<u>ATCTGAACCG</u>	<u>CCGTTCTTTG</u>	<u>CGTCTTGATG</u>
451	<u>TTAGCCC</u>				

Key

single bold underline	nt 1~121 = goose amplicon (121 bp)
double underline	nt 122~226 = duck amplicon (105 bp)
dashed underline	nt 227~344 = turkey amplicon (118 bp)
dotted underline	nt 345~457 = pigeon amplicon (113 bp)

Figure 2 — Complete sequence of nucleotides (nt) and annotation of the insertion in plasmid pUC57

Each participant in the collaborative trial received a solution containing plasmid pUC57 DNA adjusted to 1 000 copies/ μ l of the target sequence (see [Figures 1](#) and [2](#)) in 20 ng/ μ l sonicated salmon sperm DNA. The concentration was measured before distribution by digital PCR (QX100 Droplet Digital PCR System²⁾). Serial dilutions were produced by the 12 laboratories in the range of 0,02 copies/ μ l to 4 copies/ μ l using 0,2 $\times$ TE buffer containing 20 ng/ μ l sonicated salmon sperm DNA. Each participant measured six replicates per concentration level. A positive result was achieved for five copies per PCR in 72 out of 72 tests (see [Table 5](#)).

Probability of detection (POD) describes the probability that PCR amplification will take place at a given number of copies of the target sequences (see ISO/TS 16393^[4]). Qualitative data generated across all laboratories and dilution levels (see [Table 5](#)) was used to determine the POD = 0,95 of the detection method (see [Table 6](#)) as described in Reference [\[3\]](#). Standard deviation was determined to be 0,30 and the LOD_{95 %} was 3,1 copies; both parameters well below the required maximum of 1 and 20 copies, respectively^[5].

Table 5 — Collaborative trial results for the limit of detection (LOD_{95 %})

Copy number of specific DNA sequence in duck genomic sequence per PCR (nominal)	Number of positive results ($C_t < 45$) out of 72 results
20	72
10	72
5	72
2	68
1	49
0,5	34
0,1	8

2) This is a product supplied by Bio-Rad GmbH. This information is given for the convenience of users of this document and does not constitute an endorsement by ISO of the product named. Equivalent products may be used if they can be shown to lead to the same results.

Table 6 — Collaborative trial results for the probability of detection (POD)

Parameter		Specific DNA sequence in duck genomic sequence
Number of laboratories		12
Number of PCR replicates per dilution level		6
POD curve	The mean probability of detection across laboratories (LPOD)	0,80
	95 % confidence interval for LPOD	0,76 to 0,83
	Slope b relative to the ideal POD curve (b = 1)	1,16
	Laboratory standard deviation, σ_L	0,30
LOD _{95 %} (in copies)	Theoretical median laboratory	3,1

9.5 Specificity

A representative sequence from the *Anas platyrhynchos* breed pekin duck isolate CAU_Pekin_2.0 Chr1, whole genome shotgun sequence (i.e. GenBank accession number JACEUM010000001.1) was selected as a PCR target according to the alignment results among reference genomes of various bird species at the beginning of 2019.^[1] Primers and probes were designed and optimized using primer-probe selection and optimization software.

The theoretical exclusive specificity of *Anas platyrhynchos* genome sequence's primers and probes was analysed for homology to other species using the BLASTN program.^[6] The 105-bp sequence used as query is part of the NCBI accession number JACEUM010000001.1 (nucleotides position: 200 064 152 – 200 064 048). Similarity search results are given in [Annex A](#). There was no homology with other genes and species.

The assays specified in [Table 7](#) were established with DNA from different species (about 200 ng/PCR). Theoretically expected data were established by queries in public NCBI databases^[6].

This PCR amplification is specific to mallard duck (*Anas platyrhynchos*) and spot-billed duck (*Anas zonorhyncha*) DNA detection as expected. Further experimental tests observed cross detection of Mandarin duck (*Aix galericulata*), white-winged duck (*Asarcornis scutulata*), tufted duck (*Aythya fuligula*) and muscovy ducks (*Cairina moschata*).

The inclusive specificity was tested against nine breeds including *Anas platyrhynchos* Common Mallard (unknown origin, egg and meat), Khaki Campbell (England, egg), Cherry Valley (England, meat), Pekin-duck (China, meat), Tegel (Australia, meat), Shaoxing (China, egg), Jinding (China, egg), Gaoyou (China, egg and meat) and *Anas zonorhyncha* breed spot-billed (China, wild). At approximate 100 copies of target DNA, all of the breed samples were detected with the expected positive signals and amplification curves. Inclusivity of the 105 base target sequence was also evaluated using the BLASTN program against the GenBank whole animal genomes database. Results indicating that the 105 base target sequence is unique for duck are provided in [Annex A](#).

Members of the *Anatidae* family^[7] and its family tree established with available public genomic sequences are shown in Figure B.1.

Table 7 — Specificity of the target duck genomic sequence detection method

Species test		Theoretically expected	Experimental confirmation
Animal	Bison (<i>Bison bison</i>)	N	N
	Black swan (<i>Cygnus atratus</i>)	N	N
	Camel (<i>Camelus bactrianus</i>)	N	N
	Carp (<i>Cyprinus carpio</i>)	N	N
	Cat (<i>Felis catus</i>)	N	N
	Cattle (<i>Bos taurus</i>)	N	N
	Chicken (<i>Gallus gallus</i>)	N	N
	Domestic goose (<i>Anser anser domesticus</i>)	N	N
	Domestic turkey (<i>Meleagris gallopavo</i>)	N	N
	Dog (<i>Canis familiaris</i>)	N	N
	Donkey (<i>Equus asinus</i>)	N	N
	Elk (<i>Cervus canadensis</i>)	N	N
	Goat (<i>Capra hircus</i>)	N	N
	Goldfish (<i>Carassius auratus</i>)	N	N
	Horse (<i>Equus caballus</i>)	N	N
	Indian zebu (<i>Bos indicus</i>)	N	N
	Mallard duck (<i>Anas platyrhynchos</i>)	Pos	Pos
	Mandarin duck (<i>Aix galericulata</i>)	Pos	Pos
	Mouse (<i>Mus musculus</i>)	N	N
	Muscovy duck (<i>Cairina moschata</i>)	Pos	Pos
	Ostrich (<i>Struthio camelus</i>)	N	N
	Pheasant (<i>Phasianus colchicus</i>)	N	N
	Pig (<i>Sus scrofa domesticus</i>)	N	N
	Quail (<i>Coturnix coturnix</i>)	N	N
	European rabbit (<i>Oryctolagus cuniculus</i>)	N	N
	Rainbow trout (<i>Oncorhynchus mykiss</i>)	N	N
	Rat (<i>Rattus norvegicus</i>)	N	N
	Rhesus macaque (<i>Macaca mulatta</i>)	N	N
	Rock pigeon (<i>Columba livia</i>)	N	N
	Sheep (<i>Ovis aries</i>)	N	N
	Eastern spot-billed duck (<i>Anas zonorhyncha</i>)	Pos	Pos
	Swan goose (<i>Anser cygnoides domesticus</i>)	N	N
	Tufted duck (<i>Aythya fuligula</i>)	Pos	Pos ^a
	Eastern spotted dove (<i>Streptopelia chinensis</i>)	N	N
	Water buffalo (<i>Bubalus bubalis</i>)	N	N
	Wild turkey (<i>Meleagris Ocellata</i>)	N	N
	White-winged duck (<i>Asarcornis scutulata</i>)	Pos	Pos ^b

Key

Pos: positive; N: negative

^a In place of a biological sample of Tufted duck (*Aythya fuligula*) for experimental testing, the high-homology 105 base target sequence of accession number WNMM01000001.1 (199 307 524 bp – 199 307 628 bp) identified by BALSTN was synthesized and cloned in plasmid pUC57.

^b In place of a biological sample of White-winged duck (*Asarcornis scutulata*) for experimental testing, the high-homology 105 base target sequence of accession number VZS001008967.1 (752 bp – 648 bp) identified by BALSTN was synthesized and cloned in plasmid pUC57.

Table 7 (continued)

Species test		Theoretically expected	Experimental confirmation
Human	Wild yak (<i>Bos mutus</i>)	N	N
	Human (<i>Homo sapiens</i>)	N	N
Plant	Alfalfa (<i>Medicago sativa</i>)	N	N
	Corn (<i>Zea mays</i>)	N	N
	Rapeseed (<i>Brassica napus</i>)	N	N
	Rice (<i>Oryza sativa</i>)	N	N
	Sorghum (<i>Sorghum bicolor</i>)	N	N
	Soya (<i>Glycine max</i>)	N	N
	Wheat (<i>Triticum aestivum</i>)	N	N

Key

Pos: positive; N: negative

^a In place of a biological sample of Tufted duck (*Aythya fuligula*) for experimental testing, the high-homology 105 base target sequence of accession number WNMM01000001.1 (199 307 524 bp – 199 307 628 bp) identified by BALSTN was synthesized and cloned in plasmid pUC57.

^b In place of a biological sample of White-winged duck (*Asarcornis scutulata*) for experimental testing, the high-homology 105 base target sequence of accession number VZS001008967.1 (752 bp – 648 bp) identified by BALSTN was synthesized and cloned in plasmid pUC57.

10 Test report

The test report should be prepared as specified in ISO 20813 and other applicable standards (e.g. ISO 24276).

Annex A

(informative)

BlastN +2.12.0 results for query of GenBank RefSeq genome (refseq_genomes) and whole-genome shotgun contigs (wgs)

A.1 Query

A.1.1 Query ID: JACEUM010000001.1 (200 064 152 to 200 064 048 bp).

A.1.2 Description: *Anas platyrhynchos* breed pekin duck isolate CAU_Pekin_2.0 Chr1, whole genome shotgun sequence.

A.1.3 Molecule type: nucleic acid.

A.1.4 Query length: 105 bp.

A.2 Descriptions

See [Table A.1](#).

Table A.1 — Descriptions

Description	Max score	Total score	Query cover%	E value	% Ident	Accession	Accession length	Standard data-bases	Organism
<i>Anas platyrhynchos</i> breed pekin duck isolate CAU_Pekin_2.0 Chr1, whole genome shotgun sequence	195	195	100	2e-47	100	JACE-UM010000001	207246783	Whole-genome Shotgun Contigs (wgs)	birds (taxid: 8782)
<i>Anas zonorhyncha</i> breed spot-billed scaffold627, whole genome shotgun sequence	195	195	100	2e-47	100	NOIK01000625.1	360378	Whole-genome Shotgun Contigs (wgs)	birds (taxid: 8782)
<i>Asarcornis scutulata</i> isolate OUT-0051 scaffold113644, whole genome shotgun sequence	191	191	98	3e-46	100	VZS001008967.1	1936	Whole-genome Shotgun Contigs (wgs)	birds (taxid: 8782)
<i>Aythya fuligula</i> isolate bAytFul2 chromosome 1, bAytFul2.pri	185	185	98	9e-44	99	NC_045559.1	207018403	GenBank RefSeq Genome (refseq_genomes)	birds (taxid: 8782)
<i>Aix galericulata</i> isolate NCHU-2012 Aplatyrhynchos.1_RagTag, whole genome shotgun sequence	180	180	98	6e-43	98,06	JAL-LAJ010000001.1	201074787	Whole-genome Shotgun Contigs (wgs)	birds (taxid: 8782)
<i>Cairina moschata</i> domestica isolate CanardBarbarie3_02_2015_G Scaffold136, whole genome shotgun sequence	169	169	98	1e-39	96,12	QZEJ01000136.1	2119812	Whole-genome Shotgun Contigs (wgs)	birds (taxid: 8782)

ISO/TS 20224-10:2024(en)

Query	1	TCTTCACAAGCAGGGTCTAATGGAAGACTTGCTGGCCTGCTCTCTACTGGTGATGATGTG	60
Sbjct	752	TCTTCACAAGCAGGGTCTAATGGAAGACTTGCTGGCCTGCTCTCTACTGGTGATGATGTG	693
Query	61	GTGAGATGCGTGCTGTGCCTCTCTTCCTCTGGACCTTCTGCCAAG	105
Sbjct	692	GTGAGATGCGTGCTGTGCCTCTCTTCCTCTGGACCTTCTGCCAGG	648

Aythya fuligula isolate bAytFul2 chromosome 1, whole genome shotgun sequence

Sequence ID: WNMM01000001.1

Length: 207018403 Number of matches: 1 Range 1: 199307524 to 199307628

Score	Expect	Identities	Gaps	Strand
185 bits(100)	1e-44	102/103(99 %)	0/103(0 %)	Plus/Minus

Query	1	TCTTCACAAGCAGGGTCTAATGGAAGACTTGCTGGCCTGCTCTCTACTGGTGATGATGTG	60
Sbjct	199307628	TCTTCACAAGCAGGGTCTAATGGAAGACTTGCTGGCCTGCTCTCTACTGGTGATGATGTG	199307569
Query	61	GTGAGATGCGTGCTGTGCCTCTCTTCCTCTGGACCTTCTGCCAAG	105
Sbjct	199307568	GTGAGATGCGTGCTGTGCCTCTCTTCCTCTGGACCTTCTGCCAGG	199307524

Aix galericulata isolate NCHU-2012 Aplatyrrhynchos.1 RagTag, whole genome shotgun sequence

Sequence ID: JALLAJ010000001.1

Length: 201074787 Number of matches: 1 Range 1: 193845780 to 193845884

Score	Expect	Identities	Gaps	Strand
180 bits(97)	6e-43	101/103(98 %)	0/103(0 %)	Plus/Minus

Query	1	TCTTCACAAGCAGGGTCTAATGGAAGACTTGCTGGCCTGCTCTCTACTGGTGATGATGTG	60
Sbjct	193845884	TCTTCACAAGCAGGGTCTAATGGAAGACTTGCTGGCCTGCTCTCTGCTGGTGATGATGTG	193845825
Query	61	GTGAGATGCGTGCTGTGCCTCTCTTCCTCTGGACCTTCTGCCAAG	105
Sbjct	193845824	GTGAGATGCGTGCCGTGCCTCTCTTCCTCTGGACCTTCTGCCAGG	193845780

Cairina moschata domestica isolate Canard Barbarie3_02_2015_G Scaffold136, whole genome shotgun sequence

Sequence ID: QZEJ01000136.1

ISO/TS 20224-10:2024(en)

Length: 2119812	Number of matches: 1		Range 1: 1872455 to 1872559	
<u>Score</u>	<u>Expect</u>	<u>Identities</u>	<u>Gaps</u>	<u>Strand</u>
169 bits(91)	1e-39	99/103(96 %)	0/103(0 %)	Plus/Minus

Query	1	TCTTCACAAGCAGGGTCTAATGGAAGACTTGCTGGCCTGCTCTCTACTGGTGATGATGTG	60
Sbjct	1872559	TCTTCACAAGCAGGGTCTAATGGAAGACTTGCTGGCCTGCTCTCTGCTGCTGATGACGTG	1872500
Query	61	GTGAGATGCGTGCTGTGCCTCTCTTCCTCTGGACCTTCTGCCAAG	105
Sbjct	1872499	GTGAGATGCGTGCCGTGCCTCTCTTCCTCTGGACCTTCTGCCAGG	1872455

STANDARDSISO.COM : Click to view the full PDF of ISO/TS 20224-10:2024

Annex B

(informative)

Members of the *Anatidae* family and its family tree established with available public genomic sequences

B.1 Members of the *Anatidae* family

Table B.1 — Member of the *Anatidae* family

Genus No.	Genus	Species	Accession number of available public genomic sequence
1	<i>Aix</i>	<i>galericulata</i>	JALLAJ0000000000.1 (wgs ^a)
		<i>sponsa</i>	
2	<i>Alopochen</i>	<i>aegyptiaca</i>	
		<i>kervazoi</i>	
		<i>mauritiana</i>	
3	<i>Amazonetta</i>	<i>brasiliensis</i>	
^a wgs: whole genome shotgun			

Table B.1 (continued)

Genus No.	Genus	Species	Accession number of available public genomic sequence
4	Anas	<i>acuta</i>	
		<i>albogularis</i>	
		<i>andium</i>	
		<i>aucklandica</i>	
		<i>bahamensis</i>	
		<i>bernieri</i>	
		<i>capensis</i>	
		<i>carolinensis</i>	
		<i>castanea</i>	
		<i>chlorotis</i>	
		<i>crecca</i>	
		<i>diazi</i>	
		<i>eatoni</i>	
		<i>erythrorhyncha</i>	
		<i>flavirostris</i>	
		<i>fulvigula</i>	
		<i>georgica</i>	
		<i>gibberifrons</i>	
		<i>gracilis</i>	
		<i>laysanensis</i>	
		<i>luzonica</i>	
		<i>melleri</i>	
		<i>nesiotis</i>	
		<i>platyrhynchos</i>	JACEUM000000000.1 (wgs)
		<i>poecilorhyncha</i>	
		<i>rubripes</i>	
		<i>sparsa</i>	
		<i>superciliosa</i>	
		<i>theodori</i>	
		<i>undulata</i>	
		<i>wyvilliana</i>	
		<i>zonorhyncha</i>	NOIK01000625.1 (wgs)
^a wgs: whole genome shotgun			

Table B.1 (continued)

Genus No.	Genus	Species	Accession number of available public genomic sequence
5	Anser	albifrons	
		anser	
		brachyrhynchus	NXHY01000202.1 (wgs)
		caerulescens	
		canagicus	
		cygnoides	NW_025927981.1 (genome)
		erythropus	
		fabalis	
		indicus	VDDG01000056.1 (wgs)
		rossii	
		serrirostris	
6	Asarcornis	scutulata	VZSO01008967.1 (wgs)
7	Aythya	affinis	
		australis	
		baeri	JAKRSJ000000000.1 ((wgs)
		collaris	
		ferina	
		fuligula	NC_045559.1 (genome)
		innotata	
		marila	JAOQIG000000000.1 (wgs)
		novaeseelandiae	
		nyroca	
		valisineria	
8	Biziura	lobata	
9	Branta	bernicla	
		canadensis	SNRU01005576.1 (wgs)
		hutchinsii	
		leucopsis	
		ruficollis	
		sandvicensis	
10	Bucephala	albeola	
		clangula	
		islandica	
11	Cairina	moschata	QZEJ01000136.1 (wgs)
12	Callonetta	leucophrys	
13	Camptorhynchus	labradorius	
14	Cereopsis	novaehollandiae	
15	Chloephaga	hybrida	
		melanoptera	
		picta	
		poliocephala	
		rubidiceps	

^a wgs: whole genome shotgun